

# DISPOSABLE BIOPSY FORCEPS

# DISPOSABLE BIOPSY FORCEPS

PIC	Part No.	Diameter of pliers cup (mm)	Opening (mm)	Effective Length (cm)	Needle or not	Coated Plastic	Channel (mm)	Handle Type	Clamp Head Classification
	SD-01-B2316-Y	2.3	6	160	NO	YES	2.8	TYPE I Rotatable	TYPE B Alligator
	SD-01-B2323-Y	2.3		230		YES	2.8		
	SD-01-B1816-Y	1.8		160		YES	2.0		
	SD-01-B2316-N	2.3		160		NO	2.8		
	SD-01-B2323-N	2.3		230		NO	2.8		
	SD-01-B1816-N	1.8		160		NO	2.0		
	SD-01-B2316J-Y *	2.3		160		YES	2.8	TYPE II non-rotatable	
	SD-01-B2323J-Y *	2.3		230		YES	2.8		
	SD-01-B1816J-Y	1.8		160		YES	2.0		
	SD-01-B2316J-N *	2.3		160		NO	2.8		
	SD-01-B2323J-N *	2.3		230		NO	2.8		
	SD-01-B1816J-N	1.8		160		NO	2.0		
	SD-01-D2316-Y	2.3	6	160	NO	YES	2.8	TYPE I Rotatable	TYPE D Oval
	SD-01-D2323-Y	2.3		230		YES	2.8		
	SD-01-D1816-Y	1.8		160		YES	2.0		
	SD-01-D2316-N	2.3		160		NO	2.8		
	SD-01-D2323-N	2.3		230		NO	2.8		
	SD-01-D1816-N	1.8		160		NO	2.0		
	SD-01-D2316J-Y *	2.3		160		YES	2.8	TYPE II non-rotatable	
	SD-01-D2323J-Y *	2.3		230		YES	2.8		
	SD-01-D1816J-Y	1.8		160		YES	2.0		
	SD-01-D2316J-N *	2.3		160		NO	2.8		
	SD-01-D2323J-N *	2.3		230		NO	2.8		
	SD-01-D1816J-N	1.8		160		NO	2.0		

■ The \* one is our top seller model

■ Other specification can be customized according to requirement

Zhejiang Soudon Medical Technology Co., Ltd

Tel : 0086-13989453632

Fax : 0086-571-89061988


Email : info@soudonmed.com

Web : www.soudonmed.com

Add : No.1 Nangong He Road, Building 4#, LinPing Zone, Hangzhou Zhejiang, 311100



Anticipate a better endoscopic experience



**Soudon Medical**  
Anticipate a better endoscopic experience

# DISPOSABLE BIOPSY FORCEPS

It is used for sampling living tissues under endoscopy, such as digestive tract and respiratory tract, or for forceps and removal of foreign bodies.

- Sterile packaging, ready to use.
- Sharp incision, moderate sampling size and high positive detection rate.
- The tong head is designed to pass through the curved Tong path. The surface is coated with super slip coating, which has little wear on the Tong track.
- The top and bottom dead points of the handle are clear and comfortable. The product line is complete and suitable for various clinical needs.

*Brightly color* coated leads clearer under endoscope.

*Mark line* — more convenient for doctors' judgement.

