

Endoscopic Multi Band Ligator

Endoscopic Multi Band Ligator

Part No.	Endoscope Outer Diameter (mm)	Materials	Number of Bands	Trigger Cord Length (mm)
SD-13-A-4	8.6~11.3	ABS+PVC+Rubber	4	1650
SD-13-A-6			6	
SD-13-A-8			8	
SD-13-A-10			10	
SD-13-A-12			12	
SD-13-B-4	9.5~13		4	
SD-13-B-6			6	
SD-13-B-8			8	
SD-13-B-10			10	
SD-13-B-12			12	
SD-13-C-4	11~14		4	
SD-13-C-6			6	
SD-13-C-8		8		
SD-13-C-10		10		
SD-13-C-12		12		

■ Other specification can be customized according to requirement

Zhejiang Soudon Medical Technology Co., Ltd

Tel : 0086-13989453632

Fax : 0086-571-89061988

Email : info@soudonmed.com

Web : www.soudonmed.com

Add : No.1 Nangong He Road, Building 4#,LinPing Zone,
Hangzhou Zhejiang,311100



Anticipate
a better endoscopic
experience



 **Soudon Medical**
Anticipate a better endoscopic experience


ENDOSCOPIC MULTI BAND LIGATOR

Used for ligating esophageal varices combined with Endoscope.

- Precise ligation, obvious release resistance, easy to operate
- Multi-bands design achieves to treat several points in one operation
- Medical rubber material with good biocompatibility and physical property.
- Bands release reminds by feeling and sound
- Square shape cross section design make ligation stronger and safer



Yellow code color ring for quick identification of the last shot



Humanized handle design make operation more easy and comfortable.