

# Disposable Endoscope Injection Needles

| PIC  | Part No.        | Catheter OD (mm) | Effective Length (cm) | Model Needle (G) | Needle Extension Length (mm) | Metal End Cap                |
|--|-----------------|------------------|-----------------------|------------------|------------------------------|------------------------------|
|   | SD-04-A2323-4 * | 2.3              | 230                   | 23               | 4                            | TYPE A Without metal end cap |
|  | SD-04-A2323-6 * | 2.3              | 230                   | 23               | 6                            | TYPE A Without metal end cap |
|  | SD-04-A2325-4   | 2.3              | 230                   | 25               | 4                            | TYPE A Without metal end cap |
|  | SD-04-A2325-6   | 2.3              | 230                   | 25               | 6                            | TYPE A Without metal end cap |
|  | SD-04-A1823-4   | 2.3              | 180                   | 23               | 4                            | TYPE A Without metal end cap |
|  | SD-04-A1823-6   | 2.3              | 180                   | 23               | 6                            | TYPE A Without metal end cap |
|  | SD-04-A1825-4   | 2.3              | 180                   | 25               | 4                            | TYPE A Without metal end cap |
|  | SD-04-A1825-6   | 2.3              | 180                   | 25               | 6                            | TYPE A Without metal end cap |
|  | SD-04-B2323-4   | 2.3              | 230                   | 23               | 4                            | TYPE B With metal end cap    |
|  | SD-04-B2323-6   | 2.3              | 230                   | 23               | 6                            | TYPE B With metal end cap    |
|  | SD-04-B2325-4   | 2.3              | 230                   | 25               | 4                            | TYPE B With metal end cap    |
|  | SD-04-B2325-6   | 2.3              | 230                   | 25               | 6                            | TYPE B With metal end cap    |
|  | SD-04-B1823-4   | 2.3              | 180                   | 23               | 4                            | TYPE B With metal end cap    |
|  | SD-04-B1823-6   | 2.3              | 180                   | 23               | 6                            | TYPE B With metal end cap    |
|  | SD-04-B1825-4   | 2.3              | 180                   | 25               | 4                            | TYPE B With metal end cap    |
|  | SD-04-B1825-6   | 2.3              | 180                   | 25               | 6                            | TYPE B With metal end cap    |

■ The \* one is our top seller model

■ Other specification can be customized according to requirement

Zhejiang Soudon Medical Technology Co., Ltd

Tel : 0086-13989453632

Fax : 0086-571-89061988

Email : info@soudonmed.com

Web : www.soudonmed.com

Add : No.1 Nangong He Road, Building 4#, LinPing Zone,  
Hangzhou Zhejiang, 311100



# Disposable Endoscope Injection Needles

Anticipate  
a better endoscopic  
experience



**Soudon Medical**  
Anticipate a better endoscopic experience


# DISPOSABLE ENDOSCOPE INJECTION NEEDLES

Used with endoscope, suitable for submucosal injection in digestive tract.

- Tip Angle 30 degrees -- sharp puncture.
- Transparent inner tube - blood return can be observed.
- Step needle -- keep the needle in the middle of the outer tube to improve the bending performance; Increase the inner diameter of the inner tube to increase the injection flow.
- Star inner tube - small push resistance and smooth push.
- One key back needle design -- one hand operation, convenient and fast.



The needle tip not rebound during injection as **stable needle seat**.



Visual blood return as **transparent outer tube**, helpful for doctors' operate and judge.

**PP material** outer tube with better support performance.  
**Careful and slow release needle.**